

SIT8008BIL3-33E

Request a Quote

SIT8008BIL3-33E Information

| | Part Number | SIT8008BIL3-33E | |
|--------------------|--------------|---|--|
| | Manufacturer | SiTIME | |
| | Category | Crystals, Oscillators, Resonators Programmable Oscillators | |
| | Description | OSC PROG LVCMOS 3.3V EN/DS 20PPM | |
| | Package | 4-SMD, No Lead | |
| | | For the pricing/inventory/lead time, please contact us | |
| For Reference Only | | Website: https://www.heisener.com E-mail: salesdept@heisener.com | |

Certified Quality

Heisener's commitment to quality has shaped our processes for sourcing, testing, shipping, and every step in between. This foundation underlies each component we sell.



SIT8008BIL3-33E Specifications

| Manufacturer Part Number | SIT8008BIL3-33E |
|-----------------------------|---------------------------------------|
| Manufacturer | SiTIME |
| Category | Crystals, Oscillators, Resonators |
| | Programmable Oscillators |
| Package | 4-SMD, No Lead |
| Series | SiT8008B |
| Туре | MEMS (Silicon) |
| Programmable Type | Programmed as Request |
| Available Frequency Range | 1MHz ~ 110MHz |
| Function | Enable/Disable |
| Output | LVCMOS |
| Voltage - Supply | 3.3V |
| Frequency Stability | - |
| Frequency Stability (Total) | ±20ppm, ±25ppm, ±50ppm |
| Operating Temperature | -40°C ~ 85°C |
| Spread Spectrum Bandwidth | - |
| Current - Supply (Max) | 4.5mA |
| Ratings | - |
| Mounting Type | Surface Mount |
| Package / Case | 4-SMD, No Lead |
| Size / Dimension | 0.197" L x 0.126" W (5.00mm x 3.20mm) |
| Height | 0.031" (0.79mm) |
| | Report errors? |
| | |

SIT8008BIL3-33E Guarantees



Quality Guarantees

We provide 90 days warranty. * If the items you received were not in perfect quality, we would be responsible for your refund or replacement, but the items must be returned in their original condition.

SERVICE EUARANTEE

Service Guarantees

We guarantee 100% customer satisfaction. Our experienced sales team and tech support team back our services to satisfy all our customers.

SIT8008BIL3-33E Payment Methods





If you have any question about SIT8008BIL3-33E, please do not hesitate to contact us! Website: https://www.heisener.com E-mail: salesdept@heisener.com